KERN PBJ 620-3M

Multifunctional laboratory balance with single-cell weighing system



Article codes, product groups & Status information	
Model article-code	PBJ 620-3M
Model-Serie	PBJ
Assortment	Scales
Segment	Laboratory balances
Product group	Precision balances
Product usage type	Main article
EAN code (Model)	4045761206847
Predecessor Model	PBJ 620-3M-2.017a
Successor Model	PBJ 600-3M
Planned start date	2018-1HJ
Customs tariff number	90160010
Model Status	Yes
CE Label	Yes
Country of Origin	PH

Technical data - Primary	
Weighing capacity	620 g
Readability	0,001 g
Repeatability	0,001 g
Linearity	± 0,002 g
Minimum piece weight	10 mg
Weighing surface (d)	112×112 mm

Technical data - Construction	
Dimension housing (W×D×H)	210×330×70 mm
Bubble Level	Yes
Material weighing plate	Stainless steel
Revolving screw feet	Yes
Casing material	Cast aluminium
Dimensions of weighing plate (W×D)	112×112 mm
Weighing surface (d)	112×112 mm
Material draft shield	Plastic

Technical data - Display	
Display height	14 mm
Backlight	Yes

Technical data - Power supply	
Mains adapter/adapter included	CH;EURO;UK;US
Input Voltage	100 V - 240 V AC, 50 / 60 Hz
Interfaces	RS-232 standard

Technical data - Measuring system	
Calibration method	Internal calibration (automatic)
Linearity	± 0,002 g
Minimum weight (USP)	2 g
Readability	0,001 g
Recommended adjusting weight	200 g (F1)
Repeatability	0,001 g
Stabilisation time under laboratory conditions	2,5 s
Warm up time	4 h

KERN PBJ 620-3M

Multifunctional laboratory balance with single-cell weighing system

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	30 °C
Minimum ambient temperature	10 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	370×580×255 mm
Gross weight	6 kg
Net weight	4,2 kg
Shipping method	Parcel service
Packaging component - by weight - paper	165 g
Packaging component - by weight - cardboard	1,517 g
Packaging component - by weight - plastics	30 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - wood	0 g
Packaging component - by weight - foam	100 g
Delivery time	1 d

Technical data - Counting data	
Counting resolution	620.000
Mininum piece weight at piece counting (Laboratory)	1 mg
Minimum piece weight	10 mg

Technical data - Verification	
Minimum load (Min)	0,1 g

Optional Services	
DAkkS Calibration	963-103

Functions

Standard



Option

