KERN 346-09

Individual weight, OIML Class M1, 50 kg, Block, stainless steel



Article codes, product groups & Status information	
Model article-code	346-09
Model-Serie	346-0
Assortment	Scales
Segment	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761117150
Planned start date	2013-1HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary	
Nominal value	50 kg
Tolerance (OIML)	2,5 g
Case aluminium	346-090-600

Technical data - Construction	
Dimensions (W×D×H)	295×158×191 mm
Material	Stainless steel

Technical data - Measuring system	
OIML class	M1
OIML conform	Yes
Density (kg/m³)	7.950
Volume (cm³)	6.944

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	295×158×191 mm
Gross weight	50 kg
Net weight	50 kg
Shipping method	Freight
Packaging component - by weight - plastics	0 g
Delivery time	30 d

Optional Services	
DAkkS Calibration	962-646

Functions

KERN 346-09

Individual weight, OIML Class M1, 50 kg, Block, stainless steel

Standard



Option

