
KERN 308-01

Milligram weights, E1, 1 mg, Platelet, stainless steel polished



Article codes, product groups & Status information	
Model article-code	308-01
Model-Serie	308
Assortment	Scales
Segment	Test weights
Product group	OIML E1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761253421
Planned start date	2026-1HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary	
Nominal value	1 mg
Tolerance (OIML)	0,003 mg
Case aluminium	317-009-400
Case plastic	347-009-400
Case wood	338-090-200

Technical data - Construction	
Material	Stainless steel polished

Technical data - Measuring system	
Nominal value	1 mg
OIML class	E1
OIML conform	Yes
Tolerance (OIML)	0,003 mg

Density (kg/m³)	8.000
Volume (cm³)	0,0001

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	3,5×3,5×3 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	499 mg
Delivery time	1 d

Optional Services	
DAkkS Calibration	962-251

Functions

Standard

KERN 308-01

Milligram weights, E1, 1 mg, Platelet, stainless steel polished



Option

