

Ergonomically optimised chair scale – with approval for professional medical use in medical diagnostics, verification optional




1


2


2


3


4

Discover more details and matching accessories online!

Features

- Verification class III (verification is optional)
- Approved as a medical device according to regulation (EU) 2017/745
- This chair scale is the ideal measuring instrument for retirement homes, rehabilitation centres and clinics with lots of patients, and can weigh overweight people
- 1 Mobile version with two steerable rollers and particularly convenient locking brakes at the rear
- With its four wheels, this chair scale provides maximum mobility in bringing the scale to the patient. This ensures more efficient use of time for the clinic staff and greater safety for the patients who can be weighed in their familiar environment
- The rollers, with their large diameter, make it easier to get over door thresholds, edges and across the gaps into elevators
- For fragile patients, the extra wide, comfortable, ergonomically-optimised seat offers safe, secure seating during weighing
- Because of the clearly contrasting black colour of the seating surface, footrests and armrests ideally suitable for dementia patients
- 2 Two foldable armrests and footrests make transfers into the chair easier. Ideal for overweight patients or barrier-free use, e.g. for transfers from a bed to the scale
- 3 Ergonomically positioned carrying handles
- Hold function: While weighing patients that are unable to sit still, a mean average weight value is determined
- BMI function to determine underweight/normal weight/surplus weight
- High resolution readability: readability [d] can be increased by one decimal place for 5 sec. by the touch of a key
- Protective working cover included with delivery
- Option: Practical mains adapter pouch to store the mains adapter which is supplied as standard. Fitting option using two Velcro fasteners, see internet

Technical data

- 3 LCD display, digit height 25 mm
- Dimensions of display device W×D×H 200×130×60 mm
- Dimensions of seating surface W×D 390×360 mm
- Overall dimensions W×D×H 1000×640×950 mm
- Rechargeable battery pack integrated, as standard, operating time up to 40 h without backlight, charging time approx. 12 h
- Battery operation possible, 6×1.5 V AA not included, operating time up to 20 h
- 4 External mains adapter as standard, with strain-relief and disconnector plug to protect the power supply components
- Net weight approx. 24 kg

*Within the EU, official verification (conformity assessment in accordance with NAWI 2014/31/EU) is mandatory by law for scales that are intended for use as a medical device. Please order this where you state the location of installation and postcode.

STANDARD



FACTORY



Model	Weighing capacity	Readability	Verification value	Mandatory by law
				Verification
	[Max]	[d]	[e]	M III
KERN	kg	kg	kg	KERN
MCC 250K100M	250	0,1	0,1	965-129