KERN MBF 100K-2

Multifunctional baby and toddler scales with parent-child function



Article codes, product groups & Status information	
Model article-code	MBF 100K-2
Model-Serie	MBF
Assortment	Scales
Segment	Medical products
Product group	Baby balances
Product usage type	Main article
EAN code (Model)	4045761493827
Predecessor Model	MBD 20K-2; MBD 20K-2S5
Planned start date	2026-1HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	20 kg;100 kg
Readability	0,01 kg;0,05 kg

Technical data - Construction	
Material weighing plate	Plastic
Weighing surface (d)	600×295 mm
Dimensions of weighing plate (W×D)	600×295 mm
Dimension completely mounted (W×D×H)	600×354×91,5 mm
Material casing	Plastic

Technical data - Display	
Display height	20 mm

Technical data - Power supply	
Battery	3x 1,5V AAA
Operating time	51 h

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,15 kg;0,8 kg
Repeatability	0,04 kg;0,2 kg
Weighing units	Kg;lb
Number of weighing ranges	2
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	650×350×170 mm
Gross weight	2,6 kg
Net weight	1,6 kg
Shipping method	Parcel service
Delivery time	1 d

Optional Services	
DAkkS Calibration	963-129

Functions

Standard

