KERN MPN 200K-1HM

Integratable personal scale with height rod for wireless transfer of weighing data to EMR or EHR systems



Article codes, product groups & Status information	
Model article-code	MPN 200K-1HM
Model-Serie	MPN
Assortment	Scales
Segment	Medical products
Product group	Personal scales
Product usage type	Main article
EAN code (Model)	4045761242708
Planned start date	2022-2HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	250 kg
Readability	0,1 kg
Verification value (e)	0,1 kg

Technical data - Construction	
Material weighing plate	Steel, lacquered
Bubble Level	Yes
Revolving screw feet	Yes
Weighing surface (d)	365×360 mm
Dimensions of weighing plate (W×D)	365×360 mm
Dimension completely mounted (W×D×H)	365×570×2.134 mm
Material casing	ABS plastic
IP protection indicator	No IP protection
IP protection platform	IP65
IP protection load cell	IP65

Technical data - Display	
Display height	25 mm

Technical data - Power supply	
Battery	6×15 V AA
Operating time	20 h
Mains adapter/adapter included	EURO;UK
Charging time	8 h
Input Voltage	100 V - 240 V, 50 / 60 Hz
Operating time (Backlight off)	80 h
Rechargeable battery	Rchrg. battery optional
Interfaces	Alibi Speicher;Bluetooth (optional);KUP;RS-232 (optional);USB-Device (optional);WLAN

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,1 kg
Readability	0,1 kg
Recommended adjusting weight	250 kg (M1)
Repeatability	0,1 kg
Weighing units	Kg

KERN MPN 200K-1HM

Integratable personal scale with height rod for wireless transfer of weighing data to EMR or EHR systems

Weighing system	Strain gauge
Measuring ranges height meter	3- 80cm @ 0,1cm ± 0,5cm 89-205cm @ 0,1cm ± 0,5cm

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	10 °C
Minimum humidity	5 %


Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	1.010×290×445 mm
Gross weight	18 kg
Net weight	12 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	1,248 g
Packaging component - by weight - styrofoam	4,992 g
Delivery time	1 d

Technical data - Verification	
Verification Value E	0,1 kg
Minimum load (Min)	2 kg
Verification by KERN possible	Yes
Verification class	III
MPG class	1m

Optional Services	
DAkkS Calibration	963-129
Verification	965-129

Functions

Standard



Option



Factory

