KERN MGD 200K-1LS05

Elegant look combined with large, flat glass weighing plate to make it easier to step on and give maximum safety



Article codes, product groups & Status information	
Model article-code	MGD 200K-1LS05
Model-Serie	MGD
Assortment	Scales
Segment	Medical products
Product group	Bathroom scales
Product usage type	Main article
EAN code (Model)	4045761208377
Planned start date	2018-1HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	250 kg
Readability	0,1 kg
Dimension weighing platform (W×D×H)	380×330×27 mm

Technical data - Construction	
Weighing surface (d)	380×330 mm
Dimensions of weighing plate (W×D)	380×330 mm
Dimension completely mounted (W×D×H)	380×330×27 mm

Technical data - Display	
Display height	50 mm

Delivery time	1 d

Optional Services	
DAkkS Calibration	963-129S05

Technical data - Power supply	
Operating time	100 h

Technical data - Measuring system	
Calibration method	Calibration not possible
Linearity	± 1,35 kg
Readability	0,1 kg
Repeatability	1,35 kg
Weighing units	Kg;lb
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	-10 °C
Minimum humidity	5 %

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	410×352×215 mm
Gross weight	14 kg
Net weight	8 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	3,25 g
Packaging component - by weight - plastics	3,25 g

KERN MGD 200K-1LS05

Elegant look combined with large, flat glass weighing plate to make it easier to step on and give maximum safety

Functions

Standard

