KERN MCD 300K-1

Mobile chair scale in ergonomically-optimised design for safe, convenient weighing



Article codes, product groups & Status information	
Model article-code	MCD 300K-1
Model-Serie	MCD
Assortment	Scales
Segment	Medical products
Product group	Chair scales
Product usage type	Main article
EAN code (Model)	4045761202030
Planned start date	2018-1HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	300 kg
Readability	0,1 kg

Technical data - Construction	
Material weighing plate	Plastic
Bubble Level	Yes
Weighing surface (d)	390×360 mm
Dimensions of weighing plate (W×D)	390×360 mm
Dimension completely mounted (W×D×H)	625×990×985 mm
Material casing	Steel, lacquered

Technical data - Display	
Display height	25 mm

Technical data - Power supply	
Battery	6×15 V AA
Operating time	40 h

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,1 kg
Readability	0,1 kg
Recommended adjusting weight	300 kg (M1)
Repeatability	0,1 kg
Weighing units	Kg
Weighing system	Strain gauge
Max. divergency at eccentric load (1/3 max)	0,1 kg

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	10 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	1.070×690×1.080 mm
Gross weight	36 kg
Net weight	22 kg
Shipping method	Freight

Packaging component - by weight - cardboard	5,85 g
Packaging component - by weight - plastics	400 g
Packaging component - by weight - wood	7,15 g
Delivery time	2 d

Optional Services

DAkkS Calibration	963-129

Functions

Standard

