KERN 308-22

Set of milligram weights, E1, 1 mg - 500 mg, Platelet, , in plastic case



Article codes, product groups & Status information	
Model article-code	308-22
Model-Serie	T3.8-22-A
Assortment	Scales
Product Type	Test weights
Product group	OIML E1 Test weights
Product usage type	Main article
EAN code (Model)	4045761495005
Predecessor Model	308-22-2.024a
Planned start date	2026-1HJ
Customs tariff number	84239090
Country of Origin	CN
Selective distribution	0

Technical data - Primary	
Nominal value	1 mg - 500 mg
Weight set case	In plastic case

Technical data - Construction	
Is weight set	Yes
Number of weights in set	12

Technical data - Measuring system	
OIML class	E1
OIML conform	Yes

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	136×53×22 mm
Gross weight	0,10 kg
Net weight	0,10 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	0 g
Delivery time	1 d

Optional Services	
DAkkS Calibration	962-250

Functions

KERN 308-22

Set of milligram weights, E1, 1 mg - 500 mg, Platelet, , in plastic case

Standard

1 DAY

Option

DAkkS