KERN CM 60-2N

Pocket balance with integrated pocket calculator



Article codes, product groups & Status information	
Model article-code	CM 60-2N
Model-Serie	CM
Assortment	Scales
Product Type	Laboratory balances
Product group	Pocket balances
Product usage type	Main article
EAN code (Model)	4045761027381
Planned start date	2004-1HJ
Customs tariff number	90160010
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	60 g
Readability	0,01 g

Technical data - Construction	
Casing material	Plastic
Dimension housing (W×D×H)	85×130×25 mm
Material weighing plate	Plastic
Weighing surface (d)	70×80 mm

Technical data - Display	
Display height	12 mm
Backlight	Yes

Technical data - Power supply	
Battery	2×15 V AAA
Operating time	33 h

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,03 g
Weighing capacity	60 g
Readability	0,01 g
Recommended adjusting weight	50 g (M2)
Repeatability	0,02 g
Stabilisation time under laboratory conditions	2 s
Tare range	60 g
Warm up time	30 min
Weighing system	Strain gauge
Weighing units	G;gn;tl (Tw);tol

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	135×205×32 mm
Gross weight	0,35 kg
Net weight	0,20 kg

Shipping method	Parcel service
Packaging component - by weight - cardboard	454 g
Packaging component - by weight - styrofoam	34 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-127

Functions

Standard



Option

