KERN MFB 150K100

Design body analysis scale to determine body weight, body fat, body fluid, muscle mass, proportion of bone, BMI and metabolic rate/metabolic rate when exercising



Article codes, product groups & Status information	
Model article-code	MFB 150K100
Model-Serie	MGE
Assortment	KERN
Product Type	Medical products
Product group	Bathroom scales
Product usage type	Spare part;Trading parts list element;Main item
EAN code (Model)	4045761202108
Predecessor Model	MFB 150K100-2.017a
Planned start date	2018-1HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	182 kg
Readability	0,1 kg
Dimension weighing platform (W×D×H)	300×300×25 mm

Technical data - Construction	
Material weighing plate	Glass
Weighing surface (d)	300×300 mm
Dimensions of weighing plate (W×D)	300×300 mm
Dimension completely mounted (W×D×H)	300×300×25 mm

Technical data - Display	
Display height	35 mm

Technical data - Power supply	
Battery	2xCR2.032
Operating time	100 h

Technical data - Measuring system	
Calibration method	Calibration not possible
Linearity	± 0,5 kg
Readability	0,1 kg
Repeatability	0,1 kg
Weighing units	Kg;lb
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	-10 °C
Minimum humidity	5 %

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	330×330×60 mm
Gross weight	2,0 kg
Net weight	1,8 kg

Shipping method	Parcel service
UN number	3091
Packaging component - by weight - paper	39 g
Packaging component - by weight - cardboard	450 g
Packaging component - by weight - plastics	7 g
Delivery time	1 d

Optional Services	
DAkkS Calibration	963-129

Functions

Standard

