KERN MPI 200K-1

Personal scale with particularly large, flat weighing plate to make it easier to step on and give maximum safety



Article codes, product groups & Status information	
Model article-code	MPI 200K-1
Model-Serie	MPI
Assortment	KERN
Product Type	Medical products
Product group	Personal scales
Product usage type	Spare part;Trading parts list element;Main item
EAN code (Model)	4045761203310
Planned start date	2018-1HJ
Customs tariff number	84231090
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	200 kg
Readability	0,1 kg

Technical data - Construction	
Material weighing plate	Plastic
Weighing surface (d)	430×370 mm
Dimensions of weighing plate (W×D)	430×370 mm
Dimension completely mounted (W×D×H)	434×370×40 mm
Material casing	Plastic

Technical data - Display	
Display height	30 mm

Technical data - Power supply	
Battery	4x 1,5V AAA
Operating time	20 h

Technical data - Measuring system	
Calibration method	Calibration not possible
Linearity	± 1,1 kg
Readability	0,1 kg
Repeatability	0,3 kg
Weighing units	Kg;lb
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	-10 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	510×420×83 mm
Gross weight	3,6 kg
Net weight	3,2 kg
Shipping method	Parcel service
Packaging component - by weight - paper	27 g
Packaging component - by weight - cardboard	826 g

KERN MPI 200K-1

Personal scale with particularly large, flat weighing plate to make it easier to step on and give maximum safety

Packaging component - by weight - plastics	16 g
Packaging component - by weight - styrofoam	46 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-129

Functions

Standard

