KERN 326-06

Individual weight, OIML Class F1, 50 g, Eco Shape, stainless steel polished



Article codes, product groups & Status information	
Model article-code	326-06
Model-Serie	326
Assortment	Scales
Segment	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761247277
Predecessor Model	326-6; 326-6-2.018a; 327-66; 329-6
Planned start date	2026-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary	
Nominal value	50 g
Tolerance (OIML)	0,3 mg
Case aluminium	317-060-600
Case plastic	347-070-400
Case wood	317-060-100

Technical data - Construction	
Material	Stainless steel polished

Technical data - Measuring system	
Nominal value	50 g
OIML class	F1
OIML conform	Yes
Tolerance (OIML)	0,3 mg
Density (kg/m³)	7.960
Volume (cm³)	6,25

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	20×20×30 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services	
DAkkS Calibration	962-436

Functions

KERN 326-06

Individual weight, OIML Class F1, 50 g, Eco Shape, stainless steel polished

Standard



Option

