KERN MBD 20K-2

Ergonomically-optimised baby scale with large lying surface



Article codes, product groups & Status information	
Model article-code	MBD 20K-2
Model-Serie	MBD
Assortment	KERN
Product Type	Medical products
Product group	Baby balances
Product usage type	Spare part;Trading parts list element;Main item
EAN code (Model)	4045761201972
Predecessor Model	MBE 20K10
Successor Model	MBF 100K-2
Planned start date	2018-1HJ
Customs tariff number	84231090
Model Status	Yes
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	15 kg
Readability	0,01 kg

Technical data - Construction	
Material weighing plate	Plastic
Weighing surface (d)	560×225 mm
Dimensions of weighing plate (W×D)	560×225 mm
Dimension completely mounted (W×D×H)	597×350×74 mm
Material casing	Plastic

Technical data - Display	
Display height	30 mm

Technical data - Power supply	
Battery	2×15 V AAA
Operating time	20 h

Technical data - Measuring system	
Calibration method	Calibration not possible
Linearity	± 0,05 kg
Repeatability	0,01 kg
Weighing units	Kg;lb
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	80 %
Maximum operating temperature	40 °C
Minimum ambient temperature	10 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	654×480×177 mm
Gross weight	6 kg
Net weight	2,0 kg
Shipping method	Parcel service

Packaging component - by weight - paper	36 g
Packaging component - by weight - cardboard	850 g
Packaging component - by weight - plastics	12 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-128

Functions

Standard

