KERN 326-05

Individual weight, OIML Class F1, 20 g, Eco Shape, stainless steel polished



Article codes, product groups & Status information	
Model article-code	326-05
Model-Serie	326
Assortment	Scales
Segment	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761247260
Predecessor Model	326-5-2.018a; 326-5-2.025a; 326-5K+C; 327-65; 329-5
Planned start date	2026-2HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary	
Nominal value	20 g
Tolerance (OIML)	0,25 mg
Case aluminium	317-050-600
Case plastic	347-050-400
Case wood	317-050-100

Technical data - Construction	
Material	Stainless steel polished

Technical data - Measuring system	
Nominal value	20 g
OIML class	F1
OIML conform	Yes
Tolerance (OIML)	0,25 mg
Density (kg/m³)	7.960
Volume (cm³)	2,5

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	20×20×30 mm
Gross weight	0,05 kg
Net weight	0,05 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Optional Services	
DAkkS Calibration	962-435

Functions

Standard

KERN 326-05

Individual weight, OIML Class F1, 20 g, Eco Shape, stainless steel polished

1 DAY

Option

DAkkS