KERN 347-097-N00

Individual weight (Newton), 5 N, Slot weight (Type B), stainless steel fine turned



Article codes, product groups & Status information	
Model article-code	347-097-N00
Model-Serie	347-7-N
Assortment	Scales
Segment	Test weights
Product group	OIML M1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761366749
Planned start date	2024-2HJ
Customs tariff number	84239090
Country of Origin	DE

Technical data - Primary	
Tolerance (OIML)	25 mg

Technical data - Construction	
Material	Stainless steel fine turned

Technical data - Measuring system	
OIML class	M1
Density (kg/m³)	8.000
Volume (cm³)	59,5238

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	60×60×32 mm
Gross weight	0,50 kg
Net weight	0,50 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	6 d

Optional Services	
DAkkS Calibration	962-639N

Functions

KERN 347-097-N00

Individual weight (Newton), 5 N, Slot weight (Type B), stainless steel fine turned

Standard

6 DAYS

Option

DAkkS