KERN MBH 15K-3

Baby hanging scale for the calm and secure weighing of infants in a weighing cloth



Article codes, product groups & Status information	
Model article-code	MBH 15K-3
Model-Serie	MBH
Assortment	Scales
Segment	Medical products
Product group	Baby balances
Product usage type	Main article
EAN code (Model)	4045761503465
Planned start date	2026-2HJ
Customs tariff number	84238220
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	10 kg;15 kg
Readability	0,005 kg;0,01 kg

Technical data - Construction	
Dimension completely mounted (W×D×H)	101×25×203 mm
Material casing	Plastic
IP protection load cell	IP65

Technical data - Display	
Display height	15 mm

Technical data - Power supply	
Battery	2×15 V AAA

Delivery time	1 d

Optional Services	
DAkkS Calibration	963-128H

Operating time	100 h

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,03 kg
Recommended adjusting weight	10 kg (M1)
Repeatability	0,01 kg
Weighing units	Kg;lb
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	90 %
Maximum operating temperature	35 °C
Minimum ambient temperature	5 °C

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	225×150×70 mm
Gross weight	0,45 kg
Net weight	0,20 kg
Shipping method	Parcel service
Packaging component - by weight - cardboard	690 g
Packaging component - by weight - plastics	2 g
Packaging component - by weight - styrofoam	77 g

KERN MBH 15K-3

Baby hanging scale for the calm and secure weighing of infants in a weighing cloth

Functions

Standard

