KERN FOB 0.5K-4NS

Compact stainless steel balance with protection type IP65 against dust and water splashes



Article codes, product groups & Status information	
Model article-code	FOB 0.5K-4NS
Model-Serie	FOB-S
Assortment	Scales
Product Type	Food scales
Product group	Bench scales
Product usage type	Main article
EAN code (Model)	4045761161528
Planned start date	2015-1HJ
Customs tariff number	84238129
CE Label	Yes
Country of Origin	CN

Technical data - Primary	
Weighing capacity	500 g
Readability	0,1 g

Technical data - Construction	
Material weighing plate	Stainless steel
Weighing surface (d)	120×150 mm
Dimensions of weighing plate (W×D)	120×150 mm
Material casing	Stainless steel

Technical data - Display	
Display height	20 mm

Technical data - Power supply	
Battery	9 V Block
Operating time	20 h
Input Voltage	230 V AC, 50 Hz

Technical data - Measuring system	
Calibration method	Adjusting with external weight
Linearity	± 0,2 g
Recommended adjusting weight	500 g (M1)
Repeatability	0,1 g
Weighing units	G;lb
Number of weighing ranges	1
Weighing system	Strain gauge

Technical data - Environmental conditions	
Maximum humidity	95 %
Maximum operating temperature	35 °C
Minimum ambient temperature	10 °C
Minimum humidity	25 %

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	255×195×75 mm
Gross weight	0,90 kg
Net weight	0,80 kg
Shipping method	Parcel service
UN number	3496
Packaging component - by weight - paper	100 g

Packaging component - by weight - cardboard	63 g
Packaging component - by weight - plastics	63 g
Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - wood	0 g
Packaging component - by weight - foam	0 g
Delivery time	1 d

Optional Services

DAkkS Calibration	963-127

Functions

Standard



Option

