# KERN 347-085-N00

Individual weight (Newton), 2 N, Slot weight (Type A), stainless steel fine turned



| Article codes, product groups & Status information | |
|---|---|
| Model article-code | 347-085-N00 |
| Model-Serie | 347-5-N |
| Assortment | Scales |
| Segment | Test weights |
| Product group | OIML M1 Test weights |
| Product usage type | Main item;Accessories/Services |
| EAN code (Model) | 4045761277250 |
| Planned start date | 2022-1HJ |
| Customs tariff number | 84239090 |
| Country of Origin | DE |

| Technical data - Primary | |
|---|---|
| Tolerance (OIML) | 10 mg |
| Case plastic | 347-090-400 |

| Technical data - Construction | |
|---|---|
| Material | Stainless steel fine turned |

| Technical data - Measuring system | |
|---|---|
| OIML class | M1 |
| OIML conform | No |
| Density (kg/m$^3$) | 7.950 |
| Volume (cm$^3$) | 23,8095 |

| Technical data - Packing & Shipping | |
|---|---|
| Dimensions packaging (W×D×H) | 38×38×27 mm |
| Gross weight | 0,20 kg |
| Net weight | 0,20 kg |
| Shipping method | Parcel service |
| Packaging component - by weight - plastics | 0 g |
| Delivery time | 6 d |

| Optional Services | |
|---|---|
| DAkkS Calibration | 962-638N |

## Functions

# KERN 347-085-N00

Individual weight (Newton), 2 N, Slot weight (Type A), stainless steel fine turned

## Standard




## Option

