KERN 327-13-2027e

Individual weight, OIML Class F1, 5 kg, Knob, stainless steel polished



Article codes, product groups & Status information	
Model article-code	327-13-2027e
Model-Serie	327
Assortment	Scales
Segment	Test weights
Product group	OIML F1 Test weights
Product usage type	Main item;Accessories/Services
EAN code (Model)	4045761247185
Predecessor Model	327-13
Planned start date	2027K -1HJ
Customs tariff number	84239090
Country of Origin	CN

Technical data - Primary	
Nominal value	5 kg
Tolerance (OIML)	25 mg
Case aluminium	317-130-600
Case plastic	347-130-400
Case wood	317-130-100

Technical data - Construction	
Material	Stainless steel polished

Technical data - Measuring system	
Nominal value	5 kg
OIML class	F1
OIML conform	Yes
Tolerance (OIML)	25 mg
Density (kg/m³)	7.960
Volume (cm³)	625

Technical data - Packing & Shipping	
Dimensions packaging (W×D×H)	80×80×147 mm
Gross weight	5,0 kg
Net weight	5,0 kg
Shipping method	Parcel service
Packaging component - by weight - plastics	0 g
Delivery time	1 d

Technical data - Verification	
Verification possible	Yes

Optional Services	
DAkkS Calibration	962-443

Functions

Standard

KERN 327-13-2027e

Individual weight, OIML Class F1, 5 kg, Knob, stainless steel polished

1 DAY

Option

DAkkS